



Mediplus Capnomask™

Fast, reliable anaesthesia monitoring

The high-performance Capnomask™ offers a wide range of benefits including:

Enhanced patient safety – provides accurate monitoring to help prevent hypoxia post-surgery

Empowers fast clinical response – the unique, front dual connector features fully integrated oxygen and carbon dioxide lines, enabling fast, accurate detection of end-tidal CO₂ and supporting swift clinical response


Supports best practice – the Association of Anaesthetists of Great Britain and Ireland (AAGBI) recommends that continuous capnography should be used in all anaesthetised patients¹

Reduced risk of infection – integrated 0.22 micron filter reduces the risk of equipment cross-contamination and patient infection

Greater patient comfort – ergonomically moulded in soft, high-quality PVC for a stable, comfortable fit

Helps support faster recovery – impregnated with a pleasant, subtle vanilla scent to help reduce post-operative nausea and vomiting

Wide range of sizes – adult, paediatric and ophthalmic versions available




“The Capnomask™ has a superior attachment of the O₂ tubing to the mask and the mask quality is also superior and offers a great fit. I prefer it to any other mask on the market”

Dr Emad Aziz - Guy's and St Thomas' NHS Foundation Trust

Mediplus Capnomask™ range:

The Mediplus high performance Capnomask™ range includes the following:

- Adult mask
- Paediatric mask
- Ophthalmic mask



Ordering information

Product Code	Product Description	Filter	Box Quantity
A202Ac	Capnomask™ Adult Oxygen Mask	-	50
A202CMX	Capnomask™ Ophthalmic Oxygen Mask + CO ₂ sampling line 300cm with male/male connector + O ₂ tube 200cm	0.22µ	50
A202MX	Capnomask™ Adult Oxygen Mask + CO ₂ sampling line 300cm with male/male connector + O ₂ tube 200cm	0.22µ	50
A202PMX	Capnomask™ Paediatric Oxygen Mask + CO ₂ sampling line 300cm with male/male connector + O ₂ tube 200cm	0.22µ	50
A202CFX	Capnomask™ Ophthalmic Oxygen Mask + CO ₂ sampling line 300cm with male/female connector + O ₂ tube 200cm	0.22µ	50
A202FX	Capnomask™ Adult Oxygen Mask + CO ₂ sampling line 300cm with male/female connector + O ₂ tube 200cm	0.22µ	50
A202PFX	Capnomask™ Paediatric Oxygen Mask + CO ₂ sampling line 300cm with male/female connector + O ₂ tube 200cm	0.22µ	50
A202CMX-D	Capnomask™ DISS Connector Ophthalmic Oxygen Mask + CO ₂ sampling line 300cm with male/male connector + O ₂ tube 200cm	0.22µ	50
A202MX-D	Capnomask™ DISS Connector Adult Oxygen Mask + CO ₂ sampling line 300cm with male/male connector + O ₂ tube 200cm	0.22µ	50
A202PMX-D	Capnomask™ DISS Connector Paediatric Oxygen Mask + CO ₂ sampling line 300cm with male/male connector + O ₂ tube 200cm	0.22µ	50
A202CFX-D	Capnomask™ DISS Connector Ophthalmic Oxygen Mask + CO ₂ sampling line 300cm with male/female connector + O ₂ tube 200cm	0.22µ	50
A202FX-D	Capnomask™ DISS Connector Adult Oxygen Mask + CO ₂ sampling line 300cm with male/female connector + O ₂ tube 200cm	0.22µ	50
A202PFX-D	Capnomask™ DISS Connector Paediatric Oxygen Mask + CO ₂ sampling line 300cm with male/female connector + O ₂ tube 200cm	0.22µ	50
A220M	CO ₂ sampling line 300cm with male/male connector	-	50
A220F	CO ₂ sampling line 300cm with male/female connector	-	50

Clinical data

¹ Standards of monitoring which mandates capnography during anaesthesia, recovery and sedation, with the aim of improving patient safety. Updated statement from the AAGBI. Association of Anaesthetists of Great Britain & Ireland Safety Statement. September 2016. Association of Anaesthetists of Great Britain and Ireland. <https://www.aagbi.org/news/new-standards-monitoring-which-mandate-capnography-during-anaesthesia-recovery-and-sedation-aim>. Updated Sep 2016.

Mediplus are experts in the delivery of high performance medical devices that support positive patient outcomes in the following fields -



anaesthesia



gastroenterology



gynaecology



operating room



urology



Mediplus Ltd.
 +44 (0) 1494 551200
 help@mediplus.co.uk
 The Gateway Centre, Coronation Road,
 High Wycombe, HP12 3SU
 www.mediplus.co.uk



INVESTORS
IN PEOPLE

Accredited
Until 2019