

PUSEN

## Single Use Cysto-Nephro Videoscope

PC200-AS | PC200-AR | PC200-S | PC200-R

*Always Sterilized*  
*Always Available*  
*Always Works*

Improved Clinical Outcomes



## Straight Insertion Channel

Smooth navigation to low the risks of device damage



## Improved Surgeon Control

Manageable integrated suction button for optimal irrigation flow with superior intraoperative vision

Optical System	
Resolution	720 * 720 pixels
Field of View	90°
Depth of Field	5-80 mm
Insertion Portion	
Angle of Deflection	210° up/down
Diameter of Insertion Portion	5.0 mm
Working Length	380 mm
Minimum Instrument Channel Width	≥ 7.2 Fr
Additional Features	
Custom Buttons	PC200-AS, PC200-AR
Irrigation Tube Included	All Models



High Quality Imaging



Suction Control



High Brightness



Straight Channel

# PUSEN

Zhuhai Pusen Medical Technology Co., Ltd.

5/F Building 1, No.33, Keji San Road

High-tech Zone, Tangjiawan Town,

519085 Zhuhai, Guangdong, People's Republic of China

(86) 400-688-2023

Q Pusen

PUSEN Medical

PusenMedical